

# CHEST X-RAY

This project implements a comprehensive solution for binary classification of chest X-ray images (Normal vs. Pneumonia) using deep learning. The implementation covers data preprocessing, model development with both baseline CNN and transfer learning approaches, evaluation metrics, and explainability analysis.

## 1.ABOUT DATASET

### 1.1. Dataset Overview

- **Dataset Name :** Kaggle Chest X-Ray Images (Pneumonia) dataset
- **Total Samples:** 5,856 X-ray images
- **Classes:** 2 (Normal, Pneumonia)
- **Train/Val/Test Split:** 5,216 / 16 / 624 images
- **Class Distribution:**
  - Normal: 1,575 (26.9%)
  - Pneumonia: 4,281 (73.1%)

### 1.2. Image Characteristics

- **Resolution:** Variable (384×384 to 2048×2048 pixels)
- **Color Space:** Grayscale in medical format (converted to 3-channel RGB)
- **Quality:** Clinical-grade chest radiographs from pediatric patients
- **Projections:** Frontal views (PA and AP)
- **File Format:** JPEG (lossy compression)

### 1.3. DATA PREPROCESSING

#### Step 1: Image Normalization and resizing

- Standard Size: 224×224 pixels
- ImageNet statistics  
*Mean: [0.485, 0.456, 0.406]*  
*Std: [0.229, 0.224, 0.225]*

#### Step 2 : Data Augmentation

Augmentation	Parameter	Rationale
Horizontal Flip	$p = 0.5$	Safe for bilateral X-rays; left-right symmetry irrelevant
Rotation	$\pm 10^\circ$	Simulates patient positioning variation
Color Jitter	brightness = 0.2, contrast = 0.2	Accounts for X-ray acquisition variations

Affine Translation  $\pm 10\%$

Simulates patient movement/centering variation

---

*No augmentation for validation and test dataset.*

### **Train/Validation/Test Split Strategy**

**Strategy:** Systematic directory-based split

- Pre-existing split in Kaggle dataset (standard practice)
- Validation set balanced (50/50)
- Test set represents real-world distribution

### **Patient-Level Leakage Prevention:**

- Kaggle dataset already ensures no patient overlap between splits
- Confirmed: No duplicate images across train/val/test

## **1.4. DATALOADER CONFIGURATION:**

- Batch Size: 32 (balanced gradient stability and GPU memory)
- Shuffle: True for training, False for validation/test
- num\_workers: 0 (avoid multiprocessing issues with PIL)
- pin\_memory: True (GPU data transfer optimization)

## 2. MODEL DEVELOPMENT

- BaseLine CNN
- ResNet50

### 2.1. BaseLine CNN

#### *ARCHITECTURE DESIGN*

Input (3, 224, 224)

↓

Conv Block 1: Conv2d(3 → 32) + BatchNorm + ReLU + MaxPool

↓ (112 × 112)

Conv Block 2: Conv2d(32 → 64) + BatchNorm + ReLU + MaxPool

↓ (56 × 56)

Conv Block 3: Conv2d(64 → 128) + BatchNorm + ReLU + MaxPool

↓ (28 × 28)

Conv Block 4: Conv2d(128 → 256) + BatchNorm + ReLU + MaxPool

↓ (14 × 14)

Global Average Pooling

↓ (256)

Dropout (0.5)

↓

Linear (256 → 128) + ReLU + Dropout (0.5)

↓

Linear (128 → 2) [Output logits]

↓

Output (2 classes)

#### Model Statistics

- **Total Parameters:** 1,968,258
- **Trainable Parameters:** 1,968,258 (all unfrozen)

### 2.2 Transfer Learning with ResNet50

#### Transfer Learning Strategy

##### Phase 1: Feature Extraction (Frozen Backbone)

```
model = models.resnet50(pretrained=True)
```

```
# Freeze all layers
```

```
for param in model.parameters():
```

```
    param.requires_grad = False
```

```
# Replace final layer
```

```
model.fc = nn.Linear(2048, 2)

# Unfreeze new FC layer
for param in model.fc.parameters():
    param.requires_grad = True
```

WHY RESNET50??

Advantage	Benefit
Pretrained on ImageNet	Learns general visual features (edges, textures, shapes)
Residual Connections	Mitigates vanishing gradient; allows deeper networks
Transfer Learning Efficiency	Requires fewer samples and less training time
Proven Performance	State-of-the-art on medical imaging benchmarks

Model Statistics

- **Total Parameters:** 23,588,482
- **Trainable Parameters (Phase 1):** 2,050 (only FC layer)
- **Trainable Parameters (Phase 2):** 23,588,482 (all unfrozen)
- **Memory Footprint:** ~90 MB weights + activations

2.3 Training Configuration

Hyperparameters

Parameter	Value	Rationale
Learning Rate	0.001 (Baseline), 1e-5 (Fine-tuning)	Standard for Adam; lower LR is safer for fine-tuning
Batch Size	32	Balances gradient noise and GPU memory usage
Optimizer	Adam	Adaptive learning rate; works well in medical imaging


<b>Loss Function</b>	CrossEntropyLoss	Standard for multi-class (and binary) classification
<b>Epochs</b>	10–25	Enough for convergence while avoiding overfitting
<b>Weight Decay</b>	1e-4	L2 regularization to prevent overfitting
<b>LR Scheduler</b>	ReduceLROnPlateau	Lowens learning rate when validation loss stops improving

## 2.5 Evaluation Metrics


### Metrics Computed

#### 1. Accuracy

- Overall correctness
- Suitable when classes are balanced (validate on test set)



#### *Baseline CNN*



#### *ResNet50*

## 2. AUC-ROC (Area Under Receiver Operating Characteristic Curve)

- Probability that model ranks random positive higher than random negative
- Robust to class imbalance
- Primary metric for model selection

## 3. Sensitivity (Recall / True Positive Rate)

- Proportion of actual positives correctly identified
- **Critical for medical diagnosis:** Minimizes false negatives (missed pneumonia cases)
- Target:  $\geq 95\%$  (in clinical practice)

## 4. Specificity (True Negative Rate)

- Proportion of actual negatives correctly identified
- Minimizes unnecessary interventions (false alarms)
- Target:  $\geq 90\%$

## 5. Confusion Matrix

- Visual representation of classification performance
- Reveals specific error patterns

